

Product Datasheet

Anti-BOVINE IgG (H&L) (GOAT) Antibody Peroxidase Conjugated (2mg) [Out of stock] (orb2190966)

Description	Anti-BOVINE IgG (H&L) (GOAT) Antibody Peroxidase Conjugated (2mg) [Out of stock]
Species/Host	Goat
Conjugation	HRP
Tested Applications	ELISA, IHC, WB
Immunogen	Bovine IgG whole molecule
Preservatives	Lyophilized
Form/Appearance	Lyophilized
Storage	2°C to 8°C -20°C or -80°C
Note	For research use only
Application notes	Application Info: ELISA: 1:10,000 - 1:200,000 WB: 1:5,000 - 1:100,000 IHC: 1:500 - 1:5,000
Clonality	Polyclonal
MW	40000
Expiration Date	12 months from date of receipt.