

## Product Datasheet

# Anti-BOVINE IgG (H&L) (GOAT) Antibody Peroxidase Conjugated (2mg) [Out of stock] (orb2190967)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Anti-BOVINE IgG (H&L) (GOAT)<br>Antibody Peroxidase Conjugated<br>(2mg) [Out of stock]           |
| <b>Species/Host</b>        | Goat   |
| <b>Conjugation</b>         | HRP  |
| <b>Tested Applications</b> | ELISA, IHC, WB   |
| <b>Immunogen</b>           | Bovine IgG whole molecule  |
| <b>Preservatives</b>       | Lyophilized  |
| <b>Form/Appearance</b>     | Lyophilized  |
| <b>Storage</b>             | 2°C to 8°C -20°C or -80°C  |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Application Info: ELISA: 1:10,000<br>- 1:200,000 WB: 1:5,000 -<br>1:100,000 IHC: 1:500 - 1:5,000 |
| <b>Clonality</b>           | Polyclonal   |
| <b>MW</b>                  | 40000  |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |