

Product Datasheet

Human IgG/IgA/IgM (H&L) Antibody - FITC Conjugated [Out of stock] (orb2190991)

Description	Human IgG/IgA/IgM (H&L) Antibody - FITC Conjugated [Out of stock]
Species/Host	Goat
Conjugation	FITC
Tested Applications	DOT, FC, FLISA, IF, WB
Immunogen	Human IgG, IgA and IgM whole molecules
Preservatives	Lyophilized
Form/Appearance	Lyophilized
Storage	2°C to 8°C -20°C or -80°C
Note	For research use only
Application notes	Application Info: FLISA: 1:10,000 - 1:50,000 IF: 1:1,000 - 1:5,000 OTHER: FLOW CYTOMETRY 1:500 - 1:2,500
Clonality	Polyclonal
MW	390
Expiration Date	12 months from date of receipt.