

Product Datasheet

Anti-BOVINE IgG F(c) (GOAT) Antibody Fluorescein Conjugated (2mg) [Out of stock] (orb2191138)

Description	Anti-BOVINE IgG F(c) (GOAT) Antibody Fluorescein Conjugated (2mg) [Out of stock]
Species/Host	Goat
Conjugation	FITC
Tested Applications	FC, FLISA, IF
Immunogen	Bovine IgG F(c) fragment
Preservatives	Lyophilized
Form/Appearance	Lyophilized
Storage	2°C to 8°C -20°C or -80°C
Note	For research use only
Application notes	Application Info: FLISA: 1:10,000 - 1:50,000 IF: 1:1,000 - 1:5,000 OTHER: FLOW CYTOMETRY 1:500 - 1:2,500
Clonality	Polyclonal
MW	390
Expiration Date	12 months from date of receipt.