

Product Datasheet

Rat IgG F(ab')₂ Antibody - AP Conjugated [Out of stock] (orb2191203)

Description	Rat IgG F(ab') ₂ Antibody - AP Conjugated [Out of stock]
Conjugation	AP
Tested Applications	ELISA, IHC, WB
Immunogen	Rat IgG F(ab') ₂ fragment
Preservatives	Liquid (sterile filtered)
Form/Appearance	Liquid (sterile filtered)
Storage	2°C to 8°C
Note	For research use only
Application notes	Application Info: ELISA: 1:2,000 - 1:10,000 WB: 1:500 - 1:2,500 IHC: 1:200 - 1:1,000
Clonality	Polyclonal
MW	140000
Expiration Date	12 months from date of receipt.