

## Product Datasheet

# Mouse IgG (H&L) Antibody Fab - Texas Red Conjugated (orb2191579)

<b>Description</b>	Mouse IgG (H&L) Antibody Fab - Texas Red Conjugated
<b>Conjugation</b>	Texas Red
<b>Tested Applications</b>	FC, FLISA, IF
<b>Immunogen</b>	Mouse IgG whole molecule
<b>Preservatives</b>	Lyophilized
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	2°C to 8°C -20°C or -80°C
<b>Note</b>	For research use only
<b>Application notes</b>	Application Info: FLISA: 1:10,000 - 1:50,000 IF: 1:1,000 - 1:5,000 OTHER: FLOW CYTOMETRY 1:500 - 1:2,500
<b>Clonality</b>	Polyclonal
<b>MW</b>	625
<b>Expiration Date</b>	12 months from date of receipt.