

Product Datasheet

Chicken IgG (H&L) Antibody Fab - Biotin Conjugated [Out of stock] (orb2191592)

Description	Chicken IgG (H&L) Antibody Fab - Biotin Conjugated [Out of stock]
Species/Host	Goat
Conjugation	Biotin
Tested Applications	ELISA, IHC, WB
Immunogen	Chicken IgG whole molecule
Preservatives	Lyophilized
Form/Appearance	Lyophilized
Storage	2°C to 8°C -20°C or -80°C
Note	For research use only
Application notes	Application Info: ELISA: 1:80,000-1:800,000 WB: 1:2,000 - 1:10,000 IHC: 1:1,000 - 1:5,000
Clonality	Polyclonal
MW	455
Expiration Date	12 months from date of receipt.