

## Product Datasheet

# Sheep IgG (H&L) Antibody Fab - Biotin Conjugated [Out of stock] (orb2191606)

<b>Description</b>	Sheep IgG (H&L) Antibody Fab - Biotin Conjugated [Out of stock]
<b>Species/Host</b>	Donkey
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA, IHC, WB
<b>Immunogen</b>	Sheep IgG whole molecule
<b>Preservatives</b>	Lyophilized
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	2°C to 8°C -20°C or -80°C
<b>Note</b>	For research use only
<b>Application notes</b>	Application Info: ELISA: 1:20,000 - 1:100,000 WB: 1:2,000 - 1:10,000 IHC: 1:1,000 - 1:5,000
<b>Clonality</b>	Polyclonal
<b>MW</b>	455
<b>Expiration Date</b>	12 months from date of receipt.