

Product Datasheet

Mouse IgG (H&L) Antibody Fab - FITC Conjugated [Out of stock] (orb2191662)

Description	Mouse IgG (H&L) Antibody Fab - FITC Conjugated [Out of stock]
Species/Host	Goat
Conjugation	FITC
Tested Applications	DOT, IF, Multiplex Assay, WB
Immunogen	Mouse IgG whole molecule
Preservatives	Lyophilized
Form/Appearance	Lyophilized
Storage	2°C to 8°C -20°C or -80°C
Note	For research use only
Application notes	Application Info: IF: 1:500 - 1:2,500
Clonality	Polyclonal
MW	390
Expiration Date	12 months from date of receipt.