

Product Datasheet

Sheep IgG (H&L) Antibody Fab - FITC Conjugated [Out of stock] (orb2191664)

Description	Sheep IgG (H&L) Antibody Fab - FITC Conjugated [Out of stock]
Species/Host	Donkey
Conjugation	FITC
Tested Applications	FC, FLISA, IF
Immunogen	Sheep IgG whole molecule
Preservatives	Lyophilized
Form/Appearance	Lyophilized
Storage	2°C to 8°C -20°C or -80°C
Note	For research use only
Application notes	Application Info: FLISA: 1:10,000 - 1:50,000 IF: 1:1,000 - 1:5,000 OTHER: FLOW CYTOMETRY 1:500 - 1:2,500
Clonality	Polyclonal
MW	390
Expiration Date	12 months from date of receipt.