

Product Datasheet

Human IgG (H&L) Antibody Fab - FITC Conjugated [Out of stock] (orb2191669)

Description	Human IgG (H&L) Antibody Fab - FITC Conjugated [Out of stock]
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	FC, FLISA, IF
Immunogen	Human IgG whole molecule
Preservatives	Lyophilized
Form/Appearance	Lyophilized
Storage	2°C to 8°C -20°C or -80°C
Note	For research use only
Application notes	Application Info: FLISA: 1:10000~50000IF: 1:1000~5000
Isotype	IgG Fab
Clonality	Polyclonal
MW	390
Expiration Date	12 months from date of receipt.