

## Product Datasheet

# Bovine IgG (H&L) Antibody F(ab')<sub>2</sub> - FITC Conjugated [Out of stock] (orb2191754)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Bovine IgG (H&L) Antibody F(ab') <sub>2</sub> - FITC Conjugated [Out of stock]                              |
| <b>Conjugation</b>         | FITC  |
| <b>Tested Applications</b> | FC, FLISA, IF   |
| <b>Immunogen</b>           | Bovine IgG whole molecule   |
| <b>Preservatives</b>       | Lyophilized   |
| <b>Form/Appearance</b>     | Lyophilized   |
| <b>Storage</b>             | 2°C to 8°C -20°C or -80°C   |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Application Info: FLISA:<br>1:10,000 - 1:50,000 IF: 1:1,000 - 1:5,000 OTHER: FLOW CYTOMETRY 1:500 - 1:2,500 |
| <b>Clonality</b>           | Polyclonal  |
| <b>MW</b>                  | 390   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |