

Product Datasheet

Rat IgG (H&L) Antibody F(ab')₂ - Texas Red Conjugated (orb2191763)

Description	Rat IgG (H&L) Antibody F(ab') ₂ - Texas Red Conjugated
Species/Host	Goat
Conjugation	Texas Red
Tested Applications	IF, WB
Immunogen	Rat IgG whole molecule
Preservatives	Lyophilized
Form/Appearance	Lyophilized
Storage	2°C to 8°C -20°C or -80°C
Note	For research use only
Application notes	Application Info: IF: 1:500-1:2,500 OTHER: Flow Cytometry: 1:2,000-1:10,000
Clonality	Polyclonal
MW	625
Expiration Date	12 months from date of receipt.