

## Product Datasheet

# Dog IgG (H&L) Antibody F(ab')<sub>2</sub> - Texas Red Conjugated [Out of stock] (orb2191770)

<b>Description</b>	Dog IgG (H&L) Antibody F(ab') <sub>2</sub> - Texas Red Conjugated [Out of stock]
<b>Conjugation</b>	Texas Red
<b>Tested Applications</b>	FC, FLISA, IF, WB
<b>Immunogen</b>	Dog IgG whole molecule
<b>Preservatives</b>	Lyophilized
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	2°C to 8°C -20°C or -80°C
<b>Note</b>	For research use only
<b>Application notes</b>	Application Info: IF: 1:500-1:2,500 OTHER: Flow Cytometry: 1:2,000-1:10,000
<b>Clonality</b>	Polyclonal
<b>MW</b>	625
<b>Expiration Date</b>	12 months from date of receipt.