

Product Datasheet

Mouse IgG/IgA/IgM (H&L) Antibody F(ab')₂ - Texas Red Conjugated [Out of stock] (orb2191774)

Description	Mouse IgG/IgA/IgM (H&L) Antibody F(ab') ₂ - Texas Red Conjugated [Out of stock]
Species/Host	Goat
Conjugation	Texas Red
Tested Applications	FC, FLISA, IF
Immunogen	Mouse IgG IgA and IgM whole molecule
Preservatives	Lyophilized
Form/Appearance	Lyophilized
Storage	2°C to 8°C -20°C or -80°C
Note	For research use only
Application notes	Application Info: FLISA: 1:10,000 - 1:50,000 IF: 1:1,000 - 1:5,000 OTHER: FLOW CYTOMETRY 1:500 - 1:2,500
Clonality	Polyclonal
MW	625
Expiration Date	12 months from date of receipt.