

Product Datasheet

Rat IgG (H&L) Antibody F(ab')₂ - Texas Red Conjugated (orb2191777)

Description	Rat IgG (H&L) Antibody F(ab') ₂ - Texas Red Conjugated
Conjugation	Texas Red
Tested Applications	FC, FLISA, IF
Immunogen	Rat IgG whole molecule
Preservatives	Lyophilized
Form/Appearance	Lyophilized
Storage	2°C to 8°C -20°C or -80°C
Note	For research use only
Application notes	Application Info: FLISA: 1:10,000 - 1:50,000 IF: 1:1,000 - 1:5,000 OTHER: FLOW CYTOMETRY 1:500 - 1:2,500
Clonality	Polyclonal
MW	625
Expiration Date	12 months from date of receipt.