

Product Datasheet

Rabbit IgG F(c) Antibody F(ab')₂ - Biotin Conjugated (orb2191846)

Description	Rabbit IgG F(c) Antibody F(ab') ₂ - Biotin Conjugated
Species/Host	Goat
Conjugation	Biotin
Tested Applications	ELISA, IHC, WB
Immunogen	Rabbit IgG F(c) fragment
Preservatives	Lyophilized
Form/Appearance	Lyophilized
Storage	2°C to 8°C -20°C or -80°C
Note	For research use only
Application notes	Application Info: ELISA: 1:20,000 - 1:100,000 WB: 1:2,000 - 1:10,000 IHC: 1:1,000 - 1:5,000
Clonality	Polyclonal
MW	455
Expiration Date	12 months from date of receipt.