

Product Datasheet

Bovine IgG (H&L) Antibody F(ab')₂ - Biotin Conjugated [Out of stock] (orb2191885)

Description	Bovine IgG (H&L) Antibody F(ab') ₂ - Biotin Conjugated [Out of stock]
Conjugation	Biotin
Tested Applications	ELISA, IHC, WB
Immunogen	Bovine IgG whole molecule
Preservatives	Lyophilized
Form/Appearance	Lyophilized
Storage	2°C to 8°C -20°C or -80°C
Note	For research use only
Application notes	Application Info: ELISA: 1:10,000 - 1:50,000 WB: 1:1,000 - 1:5,000 IHC: 1:500 - 1:2,000
Clonality	Polyclonal
MW	455
Expiration Date	12 months from date of receipt.