

## Product Datasheet

# Guinea Pig IgG (H&L) Antibody F(ab')<sub>2</sub> - FITC Conjugated [Out of stock] (orb2192003)

<b>Description</b>	Guinea Pig IgG (H&L) Antibody F(ab') <sub>2</sub> - FITC Conjugated [Out of stock]
<b>Species/Host</b>	Goat
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	FC, FLISA, IF
<b>Immunogen</b>	Guinea Pig IgG whole molecule
<b>Preservatives</b>	Lyophilized
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	2°C to 8°C -20°C or -80°C
<b>Note</b>	For research use only
<b>Application notes</b>	Application Info: FLISA: 1:10,000 - 1:50,000 IF: 1:1,000 - 1:5,000 OTHER: FLOW CYTOMETRY 1:500 - 1:2,500
<b>Clonality</b>	Polyclonal
<b>MW</b>	390
<b>Expiration Date</b>	12 months from date of receipt.