

Product Datasheet

Rat IgG (H&L) Antibody F(ab')₂ - AP Conjugated (orb2192026)

Description	Rat IgG (H&L) Antibody F(ab') ₂ - AP Conjugated
Species/Host	Donkey
Conjugation	AP
Tested Applications	ELISA, IHC, WB
Immunogen	Rat IgG whole molecule
Preservatives	Liquid (sterile filtered)
Form/Appearance	Liquid (sterile filtered)
Storage	2°C to 8°C
Note	For research use only
Application notes	Application Info: ELISA: 1:2,000-1:10,000 WB: 1:500-1:2,500 IHC: 1:200-1:1,000
Clonality	Polyclonal
MW	140000
Expiration Date	12 months from date of receipt.