

## Product Datasheet

# Rat IgG (H&L) Antibody F(ab')<sub>2</sub> - AP Conjugated (orb2192038)

<b>Description</b>	Rat IgG (H&L) Antibody F(ab') <sub>2</sub> - AP Conjugated
<b>Species/Host</b>	Goat
<b>Conjugation</b>	AP
<b>Tested Applications</b>	ELISA, IHC, WB
<b>Immunogen</b>	Rat IgG whole molecule
<b>Preservatives</b>	Liquid (sterile filtered)
<b>Form/Appearance</b>	Liquid (sterile filtered)
<b>Storage</b>	2°C to 8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Application Info: ELISA: 1:2,000 - 1:10,000 WB: 1:500 - 1:2,500 IHC: 1:200 - 1:1,000
<b>Clonality</b>	Polyclonal
<b>MW</b>	140000
<b>Expiration Date</b>	12 months from date of receipt.