


Product Datasheet

Human c-Kit ELISA Kit (orb219270)

Description	Human c-Kit ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	KIT
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P10721
Sensitivity	80 pg/ml
Entrez	3815
Expiration Date	Please enquire.



Representative
Standard
Curve