


Product Datasheet

Human CD120b ELISA Kit (orb219272)

Description	Sandwich ELISA for quantitative detection of Human CD120b concentrations.
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Target	TNFRSF1B
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Urine Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	P20333
Sensitivity	2 pg/ml
Entrez	7133
Expiration Date	Please enquire.



Representative
Standard
Curve