


Product Datasheet

Human HER2 ELISA Kit (orb219398)

Description	Human HER2 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	ERBB2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P04626
Sensitivity	10 pg/ml
Entrez	2064
Expiration Date	Please enquire.



Representative Standard Curve