


## Product Datasheet

### Human IL6 ELISA Kit (orb219453)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Human IL6 ELISA Kit  |
| <b>Reactivity</b>          | Human  |
| <b>Range</b>               | 3.9 pg/ml - 250 pg/ml  |
| <b>Tested Applications</b> | ELISA  |
| <b>Target</b>              | IL6  |
| <b>Storage</b>             | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.  |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates<br>Experiment Notes: Colorimetric |
| <b>Sample Types</b>        | Cell Culture Supernates, Serum, Plasma, Tissue Homogenates                                     |
| <b>Uniprot ID</b>          | <b>P05231</b>  |
| <b>Sensitivity</b>         | 2 pg/ml  |
| <b>Entrez</b>              | <b>3569</b>  |
| <b>Expiration Date</b>     | Please enquire.  |



Representative  
Standard  
Curve