


Product Datasheet

Human Insulin ELISA Kit (orb219458)

Description	Sandwich ELISA for quantitative detection of Human Insulin concentrations.
Reactivity	Human
Range	15.6 pmol/L - 1000 pmol/L
Tested Applications	ELISA
Target	INS
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P01308
Sensitivity	7 pmol/L
Entrez	3630
Expiration Date	Please enquire.



Representative
Standard
Curve