


Product Datasheet

Human LIF ELISA Kit (orb219469)

Description	Sandwich ELISA for quantitative detection of Human LIF concentrations.
Reactivity	Human
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	LIF
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernatates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P15018
Sensitivity	1 pg/ml
Entrez	3976
Expiration Date	Please enquire.



Representative Standard Curve