

Product Datasheet

Tobramycin Antibody [Out of stock] (orb2195091)

Description	Tobramycin Antibody [Out of stock]
Species/Host	Mouse
Conjugation	Unconjugated
Preservatives	Purified, Liquid
Form/Appearance	Purified, Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	Protein A chromatographyProduct is 0.2um filtered.
Dilution Range	Immunogen: Tobramycin BSAAffinity Constant: Not determined.
Expiration Date	12 months from date of receipt.