


Product Datasheet

Human Spasmodysin ELISA Kit (orb219539)

Description	Sandwich ELISA for quantitative detection of Human Spasmodysin concentrations.
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Target	TFF2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q03403
Sensitivity	10 pg/ml
Entrez	7032
Expiration Date	Please enquire.



Representative
Standard
Curve