

## Product Datasheet

### CCL20 Antibody [Out of stock] (orb2195595)

<b>Description</b>	CCL20 Antibody [Out of stock]
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC, IHC-P, WB
<b>Immunogen</b>	E.coli derived recombinant Human MIP-3alpha.
<b>Preservatives</b>	Lyophilized powder
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Application Info: ELISA (2 - 4 ug/ml), IHC, IHC-P (10 ug/ml), WB (0.2 - 0.4 ug/ml)
<b>Isotype</b>	IgG1,k
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">P78556</a>
<b>NCBI</b>	<a href="#">NP_004582.1</a>
<b>Expiration Date</b>	12 months from date of receipt.