


Product Datasheet

Human TWEAK ELISA Kit (orb219569)

Description	Sandwich ELISA for quantitative detection of Human TWEAK concentrations.
Reactivity	Human
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	TNFSF12
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	043508
Sensitivity	10 pg/ml
Entrez	8742, 407977
Expiration Date	Please enquire.



Representative Standard Curve