

Product Datasheet

Anti-Sheep IgG (H&L) Antibody F(ab)'2 - FITC Conjugated [Out of stock] (orb2195931)

Description	Anti-Sheep IgG (H&L) Antibody F(ab)'2 - FITC Conjugated [Out of stock]
Species/Host	Donkey
Conjugation	FITC
Tested Applications	FC, IEM
Immunogen	Purified Sheep IgG, whole molecule
Preservatives	Clear, fluorescent yellow liquid
Form/Appearance	Clear, fluorescent yellow liquid
Storage	2°C to 8°C
Note	For research use only
Application notes	Application Info: This conjugate is suitable for immunomicroscopy and flow cytometry. The optimal working dilution should be determined by the investigator. Suggested starting dilution(s): 1:20 - 1:2,000 for most applications
Clonality	Polyclonal
Purity	> 90% based on SDS-PAGE Small amounts of intact IgG may be present.
Expiration Date	12 months from date of receipt.