

## Product Datasheet

# Anti-Sheep IgG (H&L) Antibody F(ab)'2 - Biotin Conjugated (orb2195932)

<b>Description</b>	Anti-Sheep IgG (H&L) Antibody F(ab)'2 - Biotin Conjugated
<b>Species/Host</b>	Donkey
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	FC, IEM
<b>Immunogen</b>	Purified Sheep IgG, whole molecule
<b>Preservatives</b>	Clear, colorless liquid
<b>Form/Appearance</b>	Clear, colorless liquid
<b>Storage</b>	2°C to 8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Application Info: This conjugate is suitable for immunomicroscopy and flow cytometry. The optimal working dilution should be determined by the investigator. Suggested starting dilution(s): 1:20 - 1:2,000 for most applications
<b>Clonality</b>	Polyclonal
<b>Purity</b>	> 90% based on SDS-PAGE Small amounts of intact IgG may be present.
<b>Expiration Date</b>	12 months from date of receipt.