

Product Datasheet

Anti-Rat IgG (H&L) Antibody F(ab)'2 (Cross Adsorbed) (orb2195939)

Description	Anti-Rat IgG (H&L) Antibody F(ab)'2 (Cross Adsorbed)
Species/Host	Donkey
Conjugation	HRP
Tested Applications	Array, ChIP, DOT, FC, IHC, IP, MS
Immunogen	Purified Rat IgG, whole molecule
Preservatives	Lyophilized
Form/Appearance	Lyophilized
Storage	2°C to 8°C
Note	For research use only
Application notes	Application Info: This antibody is suitable for all immunoassay applications. The optimal working dilution should be determined by the investigator.
Clonality	Polyclonal
Purity	> 90% based on SDS-PAGE Small amounts of intact IgG may be present.
Expiration Date	12 months from date of receipt.