
Product Datasheet

Anti-Rat IgG (H&L) Antibody F(ab)'₂ (Cross Adsorbed) (orb2195942)

Description

Anti-Rat IgG (H&L) Antibody F(ab)'2

Species/Host	Donkey
Conjugation	Unconjugated
Tested Applications	FC, IEM
Preservatives	Clear, colorless liquid. 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 0.05% (w/v) Sodium Azide.
Form/Appearance	Clear, colorless liquid. 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 0.05% (w/v) Sodium Azide.
Storage	Storage: 2-8C Shelf Life: Expires 1 year from date of receipt. Prepare working dilution prior to use and then discard.
Note	For research use only
Application notes	Application Info: This conjugate is suitable for immunomicroscopy and flow cytometry. The optimal working dilution should be determined by the investigator. Suggested starting dilution(s): 1:20 - 1:2,000 for most applications
Clonality	Polyclonal
Purity	> 90% based on SDS-PAGE Small amounts of intact IgG may be present.
Expiration Date	12 months from date of receipt.