

Product Datasheet

Anti-Rat IgG (H&L) Antibody F(ab)'2 - TRITC Conjugated (orb2195944)

Description	Anti-Rat IgG (H&L) Antibody F(ab)'2 - TRITC Conjugated
Species/Host	Goat
Conjugation	TRITC
Tested Applications	FC, IEM
Immunogen	Purified Rat IgG, whole molecule
Preservatives	Lyophilized
Form/Appearance	Lyophilized
Storage	2°C to 8°C
Note	For research use only
Application notes	Application Info: This conjugate is suitable for immunomicroscopy and flow cytometry. The optimal working dilution should be determined by the investigator. Suggested starting dilution(s): 1:20 - 1:2,000 for most applications
Clonality	Polyclonal
Purity	> 90% based on SDS-PAGE Small amounts of intact IgG may be present.
Expiration Date	12 months from date of receipt.