
Product Datasheet

Anti-Rat IgG (H&L) Antibody F(ab)'2 - FITC Conjugated (orb2195945)

Description

Anti-Rat IgG (H&L) Antibody F(ab)'2

Species/Host

Goat

Conjugation

FITC

Tested

FC, IEM

Applications**Immunogen**

Purified Rat IgG, whole molecule

Preservatives

Clear, fluorescent yellow liquid

Form/Appearance

Clear, fluorescent yellow liquid

Storage

Storage: Store undiluted liquid at 2-8C. For storage at -2°C, dilute with an equal volume of glycerol to prevent loss of enzymatic activity. Prepare working dilutions prior to use and then discard. Shelf Life: Expires 1 year from date of receipt. Prepare working dilution prior to use and then discard.

Note

For research use only

Application notes

Application Info: This conjugate is suitable for immunomicroscopy and flow cytometry. The optimal working dilution should be determined by the investigator. Suggested starting dilution(s): 1:20 - 1:2,000 for most applications

Clonality

Polyclonal

Purity

> 90% based on SDS-PAGE Small amounts of intact IgG may be present.

Expiration Date

12 months from date of receipt.