
Product Datasheet

Anti-Rat IgG (H&L) Antibody F(ab)'2 - Biotin Conjugated (orb2195947)

Description

Anti-Rat IgG (H&L) Antibody F(ab)'2

Species/Host	Goat
Conjugation	Biotin
Tested Applications	FC, IEM
Immunogen	Purified Rat IgG, whole molecule
Preservatives	Clear, colorless liquid
Form/Appearance	Clear, colorless liquid
Storage	Storage: Store undiluted liquid at 2-8C. For storage at -2°C, dilute with an equal volume of glycerol to prevent loss of enzymatic activity. Prepare working dilutions prior to use and then discard. Shelf Life: Expires 1 year from date of receipt. Prepare working dilution prior to use and then discard.
Note	For research use only
Application notes	Application Info: This conjugate is suitable for immunomicroscopy and flow cytometry. The optimal working dilution should be determined by the investigator. Suggested starting dilution(s): 1:20 - 1:2,000 for most applications
Clonality	Polyclonal
Purity	> 90% based on SDS-PAGE Small amounts of intact IgG may be present.
Expiration Date	12 months from date of receipt.