

## Product Datasheet

# Anti-Rat IgG (H&L) Antibody F(ab)'2 - Biotin Conjugated [Out of stock] (orb2195948)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Anti-Rat IgG (H&L) Antibody F(ab)'2 - Biotin Conjugated [Out of stock]  |
| <b>Species/Host</b>        | Goat  |
| <b>Conjugation</b>         | Biotin  |
| <b>Tested Applications</b> | FC, IEM   |
| <b>Immunogen</b>           | Purified Rat IgG, whole molecule  |
| <b>Preservatives</b>       | Clear, colorless liquid   |
| <b>Form/Appearance</b>     | Clear, colorless liquid   |
| <b>Storage</b>             | 2°C to 8°C  |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Application Info: This conjugate is suitable for immunomicroscopy and flow cytometry. The optimal working dilution should be determined by the investigator. Suggested starting dilution(s): 1:20 - 1:2,000 for most applications |
| <b>Clonality</b>           | Polyclonal  |
| <b>Purity</b>              | > 90% based on SDS-PAGE Small amounts of intact IgG may be present.   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |