


Product Datasheet

Mouse CCL20 ELISA Kit (orb219595)

Description	Mouse CCL20 ELISA Kit
Reactivity	Mouse
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Target	CCL20
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	089093
Sensitivity	1 pg/ml
Entrez	20297
Expiration Date	Please enquire.



Representative
Standard
Curve