

---

## Product Datasheet

**Anti-Rat IgG (H&L) Antibody F(ab)'2 - AP  
Conjugated [Out of stock] (orb2195950)**

**Description**

Anti-Rat IgG (H&amp;L) Antibody

stock]

**Species/Host**

Goat

**Conjugation**

AP

**Tested Applications**

FC, IEM

**Immunogen**

Purified Rat IgG, whole molecule

**Preservatives**

Clear, colorless liquid

**Form/Appearance**

Clear, colorless liquid

**Storage**

Storage: Store undiluted liquid at 2-8C. For storage at -2°C, dilute with an equal volume of glycerol to prevent loss of enzymatic activity. Prepare working dilutions prior to use and then discard. Shelf Life: Expires 1 year from date of receipt. Prepare working dilution prior to use and then discard.

**Note**

For research use only

**Application notes**

Application Info: This conjugate is suitable for immunomicroscopy and flow cytometry. The optimal working dilution should be determined by the investigator. Suggested starting dilution(s): 1:20 - 1:2,000 for most applications

**Clonality**

Polyclonal

**Purity**

> 90% based on SDS-PAGE Small amounts of intact IgG may be present.

**Expiration Date**

12 months from date of receipt.