


Product Datasheet

Mouse CD130 ELISA Kit (orb219609)

Description	Mouse CD130 ELISA Kit
Reactivity	Mouse
Range	125 pg/ml - 8000 pg/ml
Tested Applications	ELISA
Target	IL6ST
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q00560
Sensitivity	10 pg/ml
Entrez	16195
Expiration Date	Please enquire.



Representative
Standard
Curve