


Product Datasheet

Mouse CD269 ELISA Kit (orb219630)

| | |
|----------------------------|---|
| Description | Sandwich ELISA for quantitative detection of Mouse CD269 concentrations. |
| Reactivity | Mouse |
| Range | 31.2 pg/ml - 2000 pg/ml |
| Tested Applications | ELISA |
| Target | TNFRSF17 |
| Storage | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles. |
| Note | For research use only |
| Application notes | Cell Culture Supernatants, Serum, Plasma Notes: Colorimetric |
| Sample Types | Cell Culture Supernates, Serum, Plasma |
| Uniprot ID | 088472 |
| Sensitivity | 10 pg/ml |
| Entrez | 21935 |
| Expiration Date | Please enquire. |



Representative
Standard
Curve