


Product Datasheet

Mouse Inhibin beta A ELISA Kit (orb219730)

Description	Sandwich ELISA for quantitative detection of Mouse Inhibin beta A concentrations.
Reactivity	Mouse
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	INHBA
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q04998
Sensitivity	10 pg/ml
Entrez	16323
Expiration Date	Please enquire.



Representative
Standard
Curve