


Product Datasheet

Rat CCL3 ELISA Kit (orb219793)

Description	Rat CCL3 ELISA Kit
Reactivity	Rat
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	CCL3
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P50229
Sensitivity	60 pg/ml
Entrez	25542
Expiration Date	Please enquire.



Representative Standard Curve