


Product Datasheet

Rat CD178 ELISA Kit (orb219796)

Description	Sandwich ELISA for quantitative detection of Rat CD178 concentrations.
Reactivity	Rat
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	FASLG
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P36940
Sensitivity	10 pg/ml
Entrez	25385
Expiration Date	Please enquire.



Representative
Standard
Curve