

Product Datasheet

TGM2 Antibody - middle region (orb2203403)

Description	TGM2 Antibody - middle region
Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	TGM2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -2°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Application Info: Western Blot: Approx 80kDa band observed in Mouse and Rat Kidney lysates (calculated MW of 77.3kDa according to Human NP_004604.2 and 77.0kDa according to Mouse NP_033399.1). Recommended concentration: 1-3ug/ml. Peptide ELISA: antibody detection limit dilution 1:8000.
Protein Sequence	C-ERDLYLENPEIKIR
Clonality	Polyclonal
MW	77.3; 61.7
NCBI	NP_004604.2
Expiration Date	12 months from date of receipt.