
Product Datasheet

**Goat IgG Mouse IgG1 IgG2A IgG2B IgG3 IgA IgM
(HRP) antibody (orb22083)**

Description

Horseradish peroxidase-

polyclonal goat antiserum to
 Mouse IgG, IgA and Ig...

Species/Host

Goat

Conjugation

HRP

Preservatives

Peroxidase-coupled purified
 hyperimmune goat IgG lyophilized
 from a solution in phosphate
 buffered saline (PBS, pH 7.2). No
 preservative added, as it may
 interfere with the antibody
 activity.

Form/Appearance

Peroxidase-coupled purified
 hyperimmune goat IgG lyophilized
 from a solution in phosphate
 buffered saline (PBS, pH 7.2). No
 preservative added, as it may
 interfere with the antibody
 activity.

Storage

Storage: The lyophilized conjugate
 is shipped at ambient temperature
 and may be stored at +4°C;
 prolonged storage at or below -
 20°C. Spun down to remove
 insoluble particles, divided into
 small aliquots, frozen and stored
 at or below -20°C. Prior to use, an
 aliquot is thawed slowly in the
 dark at ambient temperature,
 spun down again and used to
 prepare working dilutions by
 adding sterile phosphate buffered
 saline (PBS, pH 7.2). Repeated
 thawing and freezing should be
 avoided. Working dilutions should
 be stored at +4°C, not refrozen,
 and preferably used the same day.
 If a slight precipitation occurs
 upon storage, this should be
 removed by centrifugation. It will
 not affect the performance of the
 immunoconjugate

Note

For research use only

Clonality

Polyclonal

Source

Purified polyclonal mouse IgG, and

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

 Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280
 240

Biorbyt LLC.

 68 TW Alexander Drive
Research Triangle Park
Durham, North
 Carolina
27709. United States

 Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

**Hazard
 Information**

This product is intended FOR
 RESEARCH USE ONLY and FOR