
Product Datasheet

**Goat IgG Mouse IgG1 IgG2A IgG2B IgG3 IgA IgM
(HRP) antibody (orb22086)**

Description

Horseradish peroxidase-

polyclonal goat antiserum to
Mouse IgG, IgA and Ig...

Species/Host

Goat

Conjugation

HRP

Preservatives

Peroxidase-coupled purified
hyperimmune goat IgG lyophilized
from a solution in phosphate
buffered saline (PBS, pH 7.2). No
preservative added, as it may
interfere with the antibody
activity.

Form/Appearance

Peroxidase-coupled purified
hyperimmune goat IgG lyophilized
from a solution in phosphate
buffered saline (PBS, pH 7.2). No
preservative added, as it may
interfere with the antibody
activity.

Storage

Storage: The lyophilized conjugate
is shipped at ambient temperature
and may be stored at +4°C;
prolonged storage at or below -
20°C. Spun down to remove
insoluble particles, divided into
small aliquots, frozen and stored
at or below -20°C. Prior to use, an
aliquot is thawed slowly in the
dark at ambient temperature,
spun down again and used to
prepare working dilutions by
adding sterile phosphate buffered
saline (PBS, pH 7.2). Repeated
thawing and freezing should be
avoided. Working dilutions should
be stored at +4°C, not refrozen,
and preferably used the same day.
If a slight precipitation occurs
upon storage, this should be
removed by centrifugation. It will
not affect the performance of the
immunoconjugate

Note

For research use only

Clonality

Polyclonal

Source

Purified polyclonal mouse IgG, and

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

 Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280
240

Biorbyt LLC.

 68 TW Alexander Drive
Research Triangle Park
Durham, North
Carolina
27709. United States

 Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

**Hazard
Information**

This product is intended FOR
RESEARCH USE ONLY and FOR