

---

## Product Datasheet

**Rabbit IgG Mouse IgG1 IgG2A IgG2B IgG3 IgA IgM  
(FITC) antibody (orb22099)**

**Description**

Fluorescein isothiocyanate-

polyclonal rabbit antiserum to  
Mouse IgG, IgA ...

**Species/Host**

Rabbit

**Conjugation**

FITC

**Preservatives**

FITC-coupled purified  
hyperimmune rabbit IgG  
lyophilized from a solution in  
phosphate buffered saline (PBS,  
pH 7.2). No preservative added,  
as it may interfere with the  
antibody activity.

**Form/Appearance**

FITC-coupled purified  
hyperimmune rabbit IgG  
lyophilized from a solution in  
phosphate buffered saline (PBS,  
pH 7.2). No preservative added,  
as it may interfere with the  
antibody activity.

**Storage**

Storage: The lyophilized conjugate  
is shipped at ambient  
temperature and may be stored at  
+4°C; prolonged storage at or  
below -20°C. Spun down to  
remove insoluble particles,  
divided into small aliquots, frozen  
and stored at or below -20°C. Prior  
to use, an aliquot is thawed slowly  
in the dark at ambient  
temperature, spun down again  
and used to prepare working  
dilutions by adding sterile  
phosphate buffered saline (PBS,  
pH 7.2). Repeated thawing and  
freezing should be avoided.  
Working dilutions should be stored  
at +4°C, not refrozen, and  
preferably used the same day. If a  
slight precipitation occurs upon  
storage, this should be removed  
by centrifugation. It will not affect  
the performance of the  
immunoconjugate

**Note**

For research use only

**Clonality**

Polyclonal

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

 Email: [info@biorbyt.com](mailto:info@biorbyt.com) | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280  
240

Biorbyt LLC.

 68 TW Alexander Drive<br>Research Triangle Park<br>Durham, North  
Carolina<br>27709. United States

 Email: [info@biorbyt.com](mailto:info@biorbyt.com) | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

**Hazard**

This product is intended FOR